



Manometry made simple.

Sandhill's InSIGHT Manometry System "works the way you think," setting the standard for simplicity, timeliness and adaptability in esophageal diagnostics. With InSIGHT, enjoy the benefits of state of the art esophageal manometry and gain and upgrade path to Impedance/Manometry. Featuring Sandhill's exclusive BioVIEW™ software, the system provides a choice of visual options designed to make analysis simple.

- **Guided data acquisition procedures**
- **Touch-screen user interface**
- **Application-specific displays**
- **Full-disclosure editing**

INSIGHT™ Manometry System



Guided procedures

Intuitive prompts help the user conduct manometry procedures with minimal training. The help function supports each procedure step with concise explanations. The guided procedures also ensure that the studies stay consistent, regardless of employee turnover.

Touch-screen user interface

Innovative touch-screen technology allows users to focus on the patient and quality of data collected, rather than computer equipment.

Application-specific displays

Each screen only displays input options required for that procedure step, minimizing screen clutter and confusion.

InSIGHT Roll-stand configuration includes:

- Central Unit
- Roll-stand with basket
- Touch LCD (flat panel)



Data display options

Two waveform display options include moving trace (traditional stripchart format) and moving erase bar. Waveforms can be frozen for closer inspection, or the screen split into active and frozen areas to review previously acquired waveforms while collecting new data.

Full-disclosure editing

The user can easily review and override computer-selected analysis points.

Data reporting

The InSIGHT manometry reports are highly flexible to satisfy your specific preferences. Study reports can be archived in print or digital forms. Editable report templates provide the capability to generate reports in local languages. Analyzed data can be exported to relational databases for research and data trending purposes.

BioVIEW® Software

In addition to the standard waveform view, InSIGHT's BioView software now includes ContourVIEW™, which uses contour mapping to provide a state of the art visual representation of the pressures for ease of analysis.

Systems are compatible with all solid-state probes and water infused systems.

LATITUDE®

Manometry Probes

The Cost-Effective Alternative

In addition to solid state probes and water infused systems, Sandhill offers Latitude Manometry Probes.



Latitude Esophageal and Anorectal Probes are comfortable, fast and single use.

For the first time, you can enjoy the benefits of circumferential pressure sensors in a single-use manometry probe.

Advanced Clinical Education

Comprehensive education is a critical element to the success of GI procedures and diagnostics. Under the direction of Dr. Donald Castell, world renowned esophagologist, Sandhill presents a variety of clinical training programs on GI diagnostic testing.



SANDHILL SCIENTIFIC

Sandhill Scientific, Inc
9150 Commerce Center Circle, #500
Highlands Ranch, CO 80129 USA
800.468.4556 303.470.7020
Fax: 303.470.2975

Sandhill Scientific - UK Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, WITNEY, Oxon OX29 4AN
Tel: +44 1865.882532
Fax: +44 1865.883970

Sandhill Scientific - Europe
Decopra s.r.o
Svabky 2, Budova B
180 00 Praha 8
Czech Republic
Tel: +420 266.310.111
Fax: +420 266.310.110